

HEGAR UTERINE DILATOR

double ended; Ø 19.0 / 20.0mm
 20.0cm (8")

SKU: 180-354-19
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-20
AUTHORS	HEGAR
Medicon	52.35.19
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-20
AUTHORS	HEGAR
Medicon	52.35.19

SPECIFICATIONS

[double ended](#)