

HEGAR UTERINE DILATOR

double ended; Ø 21.0 / 22.0mm
 20.0cm (8")

SKU: 180-354-21
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-22
AUTHORS	HEGAR
Medicon	52.35.21
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-22
AUTHORS	HEGAR
Medicon	52.35.21

SPECIFICATIONS

[double ended](#)