

HEGAR UTERINE DILATORS SET

with 26 pieces; Ø 1-26mm; in metal case



SKU: 180-350-99
Categories: [Open Surgery](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM192
KLS Martin	32-400-90
AUTHORS	HEGAR
Medicon	52.30.45

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM192
KLS Martin	32-400-90
AUTHORS	HEGAR
Medicon	52.30.45