

HEGAR UTERINE DILATORS SET

with 14 pieces; Ø 4-17mm; in metal case



SKU: 180-350-98
Categories: [Open Surgery](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM193
KLS Martin	32-400-84
AUTHORS	HEGAR
Medicon	52.30.44

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM193
KLS Martin	32-400-84
AUTHORS	HEGAR
Medicon	52.30.44