

## HEISS HAEMOSTATIC FORCEPS

strongly curved; serrated

20.0cm (8")




**SKU:** 040-393-20

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Bronchus Forceps</a>
<b>Aesculap</b>	BH209
<b>KLS Martin</b>	13-425-20
<b>AUTHORS</b>	<a href="#">HEISS</a>
<b>Medicon</b>	16.51.21
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Bronchus Forceps</a>
<b>Aesculap</b>	BH209
<b>KLS Martin</b>	13-425-20

## AUTHORS

[HEISS](#)

## Medicon

16.51.21

## SPECIFICATIONS

[curved](#), [serrated](#)