

HEISS HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")



SKU: 040-390-20

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH206
KLS Martin	13-422-20
AUTHORS	HEISS
Medicon	16.50.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH206
KLS Martin	13-422-20

AUTHORS

[HEISS](#)

Medicon

16.50.20

SPECIFICATIONS

[serrated](#), [straight](#)