


HENLEY SELFRETAINING RETRACTOR

blunt; 3X4 teeth; complete with blades

16.5cm (6.6")



SKU: 070-301-00
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-768-00
AUTHORS	HENLEY
SPECIFICATIONS	blunt , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-768-00
AUTHORS	HENLEY
SPECIFICATIONS	blunt , toothed