

HENLEY SELFRETAINING RETRACTOR ONLY

blunt; 3X4 teeth


16.0cm (6.4")



SKU: 070-301-16
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	HENLEY
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	HENLEY
SPECIFICATIONS	blunt , toothed