

HENLEY SELFRETAINING RETRACTOR

sharp; 3X4 teeth; complete with blades


16.5cm (6.6")



SKU: 070-300-00
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-769-00
AUTHORS	HENLEY
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-769-00
AUTHORS	HENLEY
SPECIFICATIONS	sharp , toothed