

HF BALL ELECTRODE; BIPOLAR; STERILE

Ø 2mm, single use

for 18.5 Fr. system; green; Box=5 pieces



SKU: 582-180-43B-S

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--|
| UNIT | Box |
| RANGE | Mini Resectoscopy (Midi) |
| TYPE | Electrodes |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--|
| UNIT | Box |
| RANGE | Mini Resectoscopy (Midi) |
| TYPE | Electrodes |