

HF KNIFE ELECTRODE; BIPOLAR; STERILE

angled; single use

for 16 Fr. System; black; Box = 5 pieces



SKU: 582-150-43GUB-S

Categories: [Endoscopy](#), [GUBBINI/Gyn](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Ellipse , Original
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Ellipse , Original
TYPE	Electrodes