

## HF LOOP ELECTRODE; BIPOLAR; STERILE

horizontal; single use

for 16 Fr. System; black; Box = 5 pieces



**SKU:** 582-150-44GUB-S

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">Box</a>        |
| <b>RANGE</b> | <a href="#">Original</a>   |
| <b>TYPE</b>  | <a href="#">Electrodes</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">Box</a>        |
| <b>RANGE</b> | <a href="#">Original</a>   |
| <b>TYPE</b>  | <a href="#">Electrodes</a> |