

HF LOOP ELECTRODE; BIPOLAR; STERILE

horizontal; single use

for 16 Fr. System; black; Box = 5 pieces



SKU: 582-150-44GUB-S

Categories: [Endoscopy](#), [GUBBINI/Gyn](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Original
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Original
TYPE	Electrodes