

HF LOOP ELECTRODE; BIPOLAR; STERILE

straight; single use

for 18.5 Fr. system; green; Box=5 pieces



SKU: 582-180-41B-S

Categories: [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Mini Resectoscopy (Midi)
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Mini Resectoscopy (Midi)
TYPE	Electrodes