

## HF LOOP ELECTRODE; BIPOLAR; STERILE

angled; single use; white-blue

for 19/22 Fr. System; Box = 5 pieces



**SKU:** 630-340-00B-S

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Resectoscopy 2.9</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>RANGE</b>	<a href="#">Resectoscopy 2.9</a>
<b>TYPE</b>	<a href="#">Electrodes</a>