

HF LOOP ELECTRODE; BIPOLAR; STERILE

angled; single use; white-blue

for 19/22 Fr. System; Box = 5 pieces



SKU: 630-340-00B-S

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Resectoscopy 2.9
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
RANGE	Resectoscopy 2.9
TYPE	Electrodes