

## HF-LOOP ELECTRODE

90° angled; tearform

for 16 Fr. System; bipolar



**SKU:** 582-150-40GUBE

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">Ellipse</a>    |
| <b>TYPE</b>  | <a href="#">Electrodes</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">Ellipse</a>    |
| <b>TYPE</b>  | <a href="#">Electrodes</a> |