

HIGH-FLOW HYSTEROSCOPY OPERATING SHEATH

Ø4.3mm; 5CH Instrument channel for telescopes

2mm/ WL: 302mm/ 30° with stopcocks

color: red

SKU: 640-022-53

Categories: [Endoscopy](#), [Gynecology](#)




ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Hysteroscopy](#)

TYPE [Operating Sheaths](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Hysteroscopy](#)

TYPE [Operating Sheaths](#)