

## HIGH-FLOW HYSTEROSCOPY OPERATING SHEATH

Ø4.3mm; 5CH Instrument channel for telescopes

2mm/ WL: 302mm/ 30° with stopcocks

color: red

**SKU:** 640-022-53

**Categories:** [Endoscopy](#), [Gynecology](#)




### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Hysteroscopy](#)

**TYPE** [Operating Sheaths](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Hysteroscopy](#)

**TYPE** [Operating Sheaths](#)