

HIGH-FLOW HYSTEROSCOPY OPERATING SHEATH

Ø4.3mm; 5CH instrument channel with stopcocks

for telescopes 2mm/ 260mm/ 30°

color: red


SKU: 640-022-43

Categories: [Endoscopy](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Hysteroscopy](#)
TYPE [Operating Sheaths](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Hysteroscopy](#)
TYPE [Operating Sheaths](#)