

HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 2.5 x 3000 mm Length

tube colour green

HP
HIGH PERFORMANCE LINE



SKU: 540-822-30GN

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)