


HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 2.5 x 5000 mm Length

tube colour lilac

HP
HIGH PERFORMANCE LINE



SKU: 540-822-50LL

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HP (High Performance)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HP (High Performance)