

## HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 2300 mm Length

tube colour yellow

**HP**  
HIGH PERFORMANCE LINE



**SKU:** 540-823-23YW

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Light Guide Cables](#)  
**TYPE** [HP \(High Performance\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Light Guide Cables](#)  
**TYPE** [HP \(High Performance\)](#)