

## HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 5000 mm Length

tube colour green

**HP**  
HIGH PERFORMANCE LINE



**SKU:** 540-823-50GN

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HP (High Performance)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">HP (High Performance)</a>