


HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 3.5 x 5000 mm Length

tube colour lilac

HP
HIGH PERFORMANCE LINE



SKU: 540-823-50LL

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)