


## HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 4.9 x 1800 mm Length

tube colour lilac

**HP**  
HIGH PERFORMANCE LINE



**SKU:** 540-824-18LL

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HP \(High Performance\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HP \(High Performance\)](#)