

HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 4.9 x 1800 mm Length

tube colour yellow

HP
HIGH PERFORMANCE LINE



SKU: 540-824-18YW

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#),
[Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HP \(High Performance\)](#)