

HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 4.9 x 2300 mm Length

tube colour yellow

HP
HIGH PERFORMANCE LINE



SKU: 540-824-23YW

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Light Guide Cables](#)
TYPE [HP \(High Performance\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Light Guide Cables](#)
TYPE [HP \(High Performance\)](#)