

HIGH PERFORMANCE LIGHT GUIDE CABLE

for LED light source; w/o connectors

Ø 4.9 x 2300 mm Length

tube colour yellow

HP
HIGH PERFORMANCE LINE



SKU: 540-824-23YW

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---------------------------------------|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HP (High Performance) |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---------------------------------------|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HP (High Performance) |