

HOFFMANN TUBING CLAMP

with fixation screw

17mm



SKU: 250-404-17

Categories: [Open Surgery](#), [Various Items](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-424-17
AUTHORS	HOFFMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-424-17
AUTHORS	HOFFMANN