

## HOFFMANN TUBING CLAMP

with fixation screw

25mm



**SKU:** 250-404-25

**Categories:** [Open Surgery](#), [Various Items](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>KLS Martin</b>	63-424-25
<b>AUTHORS</b>	<a href="#">HOFFMANN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>KLS Martin</b>	63-424-25
<b>AUTHORS</b>	<a href="#">HOFFMANN</a>