

HOFFMANN TUBING CLAMP

with fixation screw

25mm



SKU: 250-404-25

Categories: [Open Surgery](#), [Various Items](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-424-25
AUTHORS	HOFFMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-424-25
AUTHORS	HOFFMANN