

HOFFMANN TUBING CLAMP

with fixation screw

13mm



SKU: 250-404-13

Categories: [Open Surgery](#), [Various Items](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-424-13
AUTHORS	HOFFMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-424-13
AUTHORS	HOFFMANN