

HOHMANN BONE LEVER

8mm; narrow tip, 1 whole handle

22.0cm (8.8")


SKU: 190-572-08

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK185
KLS Martin	23-541-08
AUTHORS	HOHMANN
Medicon	77.11.10



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK185
KLS Martin	23-541-08
AUTHORS	HOHMANN

Medicon

77.11.10