

HOHMANN BONE LEVER

18mm; narrow tip, 1 whole handle

24.0cm (9.6")


SKU: 190-572-18

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK166
KLS Martin	23-541-18
AUTHORS	HOHMANN
Medicon	77.11.18



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK166
KLS Martin	23-541-18
AUTHORS	HOHMANN

Medicon

77.11.18