

## HOHMANN BONE LEVER

43mm; narrow tip

24.0cm (9.6")


**SKU:** 190-572-43

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK188
<b>KLS Martin</b>	23-543-43
<b>AUTHORS</b>	<a href="#">HOHMANN</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK188
<b>KLS Martin</b>	23-543-43
<b>AUTHORS</b>	<a href="#">HOHMANN</a>