

HOHMANN BONE LEVER

22mm; wide tip; 1 whole handle

26.0cm (10.4")


SKU: 190-573-22

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK172R
KLS Martin	23-547-24
AUTHORS	HOHMANN
Medicon	77.11.24



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK172R
KLS Martin	23-547-24
AUTHORS	HOHMANN

Medicon

77.11.24