

## HOHMANN BONE LEVER

15mm

16.0cm (6.4")

**SKU:** 190-571-15

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                             |
|-------------------|-----------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>        |
| <b>TYPE</b>       | <a href="#">Bone Levers</a> |
| <b>KLS Martin</b> | 23-539-15                   |
| <b>AUTHORS</b>    | <a href="#">HOHMANN</a>     |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                             |
|-------------------|-----------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>        |
| <b>TYPE</b>       | <a href="#">Bone Levers</a> |
| <b>KLS Martin</b> | 23-539-15                   |
| <b>AUTHORS</b>    | <a href="#">HOHMANN</a>     |