

HOOK FLEXIBLE

3 teeth; sharp

16.0cm (6.4")


SKU: 062-088-03

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
Aesculap	BT133
KLS Martin	15-078-03
Medicon	20.64.53
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
Aesculap	BT133
KLS Martin	15-078-03
Medicon	20.64.53

SPECIFICATIONS

[sharp, toothed](#)