

## HOOK FLEXIBLE

3 teeth; sharp

16.0cm (6.4")


**SKU:** 062-088-03

**Categories:** [Open Surgery](#), [Retractors](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	BT133
<b>KLS Martin</b>	15-078-03
<b>Medicon</b>	20.64.53
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	BT133
<b>KLS Martin</b>	15-078-03
<b>Medicon</b>	20.64.53

## SPECIFICATIONS

[sharp](#), [toothed](#)