

HOOK FLEXIBLE

3 teeth; blunt

16.0cm (6.4")


SKU: 062-089-03

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
Aesculap	BT137
KLS Martin	15-079-03
Medicon	20.65.53
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
Aesculap	BT137
KLS Martin	15-079-03
Medicon	20.65.53

SPECIFICATIONS

[blunt](#), [toothed](#)