


## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 2.5 x 3000 mm Length

tube colour green

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-812-30GN

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>                |
| <b>TYPE</b>  | <a href="#">HTT (High Temperature Technology)</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>                |
| <b>TYPE</b>  | <a href="#">HTT (High Temperature Technology)</a> |