

## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 3.5 x 2300 mm Length

tube colour yellow

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-813-23YW

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>                |
| <b>TYPE</b>  | <a href="#">HTT (High Temperature Technology)</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Light Guide Cables</a>                |
| <b>TYPE</b>  | <a href="#">HTT (High Temperature Technology)</a> |