


## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 2.5 x 1800 mm Length

tube colour lilac

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-812-18LL

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)