

HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 4.9 x 2300 mm Length

tube colour lilac

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-815-23LL

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HTT (High Temperature Technology) |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HTT (High Temperature Technology) |