

HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 4.9 x 3000 mm Length

tube colour yellow

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-815-30YW

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)