

HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 2.5 x 2300 mm Length

tube colour gray

HTT

HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-812-23

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HTT (High Temperature Technology) |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Light Guide Cables |
| TYPE | HTT (High Temperature Technology) |