


## HTT FIBER LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Ø 2.5 x 2300 mm Length

tube colour lilac

**HTT**  
HIGH TEMPERATURE TECHNOLOGY LINE



**SKU:** 540-812-23LL

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Guide Cables](#)

**TYPE** [HTT \(High Temperature Technology\)](#)