


HTT FIBER Y LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Y-Typ 2300mm; 3.5/4.8; 2x Proj / 1x Endo

tube colour green

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-817-23GN

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HTT (High Temperature Technology)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	HTT (High Temperature Technology)