


HTT FIBER Y LIGHT GUIDE CABLE

for Xenon light source; w/o connectors

Y-Typ 1800mm; 3.5/4.8; 1x Proj / 2x Endo

tube colour lilac

HTT
HIGH TEMPERATURE TECHNOLOGY LINE



SKU: 540-816-18LL

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#), [Arthroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Guide Cables](#)

TYPE [HTT \(High Temperature Technology\)](#)